Competence Flow[®]

Seite 2	- de:	Gebrauchsanweisung
Page 3	- en:	Instructions for use
Page 4	- fr:	Manuel d'instruction
Página 5	- es:	Instrucciones de uso
Pagina 6	- it:	Istruzioni per l'uso
Strana 7	- bs:	Uputstvo za upotrebu
Pagină 8	- ro:	Instrucțiuni de utilizare
Sayfa 9	- tr:	Kullanma kılavuzu
Strona 10	- pl:	Instrukcja użytkowania
Стр. 11	- ru:	Инструкция по применению

We know how

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Patientenzielgruppe

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Instructions for use

Patient target group

Children, adolescents, adults and older patients

Intended users

Dentist

Composition

Qualitative composition:

1,4-butandiol-dimethacrylate, Bis-GMA, aliphatic urethane dimethacrylate, butylated hydroxytoluene, camphorquinone, silicon dioxide, glass ceramic, iron oxide pigments, 2-(2H-Benzotriazol-2-yl)-p-cresol, 2-ethylhexyl-4-(di-methyl-amino) benzoate, diphenyl (2,4,6-trimethyl-benzoyl) phosphine oxide

Quantitative composition:

1,4-butandiol-dimethacrylate 5 %, aliphatic urethane dimethacrylate 22 %, Bis-GMA 18 %, silicon dioxide 6 %, activators and UV absorber <1 %, glass ceramic 48 %, pigments <1 %

Intended purpose of the product

Competence Flow is a light-curing, radiopaque, free-flowing composite filling material.

Indications

Cavities of classes III, IV and V.

Contraindications

Do not use in patients with known allergy to the ingredients.

Performance characteristics of the product

Competence Flow is a free-flowing composite that cures with light.

Specific precautions

None

Application

1. Clean and dry the preparation. Determine the colour.
2. Prepare the cavity according to the general rules of adhesive technique. Bevel the edges of the enamel.
3. Apply etching gel to the enamel edges for 20 to 35 seconds.
4. Completely rinse off the etching gel.
5. Dry with oil-free air. Avoid over-drying the dentin.
6. Apply adhesive as described in the usage instructions.
7. Select the syringe in the desired colour. Remove the syringe cap. Screw on the Luer Lock canula.
8. Competence Flow is applied directly by pressing the syringe. Apply Competence Flow in layers no thicker than 2 mm to ensure complete curing.
9. Leave each layer to light cure for 30 to 40 seconds.
10. Repeat steps 8 and 9 until the cavity is filled.
11. Remove excess, contour and finish with diamond finishers and flexible discs of decreasing grain size.

Warnings**Residual risks**

Please observe carefully

Light units should be emitting at 450 nm at 500 mW/cm² and should be checked regularly. Place the light as close as possible to the filling material. Do not remove the oxygen-inhibited layer, as this is needed for bonding with the next layer. Using a die is recommended. Follow the instructions for use of other materials used with Competence Flow.

Safety information

Use is contraindicated if drainage or the prescribed application method is not possible. Do not use directly on opened pulp or if there is a known allergy to any ingredient. No side effects have been reported to date. In individual cases, hypersensitivity to components of the material may occur. Avoid skin contact with the paste. In case of accidental contact, wash immediately with soap and water. In case of eye contact, rinse immediately with plenty of water and consult a doctor. Avoid materials containing eugenol, moisture and oily air, as they can prevent polymerisation at the point of contact. Discolouration may occur when using cationic mouthwashes, plaque disclosers or chlorhexidine.

Side effects

Adverse side effects are not known.
Risks in case of overdose are not known.

Interactions

Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage

Store in a dry place. The Product should be stored between 0° C - 22°C.

Minimum shelf life

Do not use after the expiry date.

Variants

REF WP4197 3 x syringes (A2, A3, A3,5), 1,0 ml + 9 cannulas
 REF WP4199 1 x syringe A1 with 1,0 ml + cannulas
 REF WP4190 1 x syringe A2 with 1,0 ml + cannulas
 REF WP4191 1 x syringe A3 with 1,0 ml + cannulas
 REF WP4192 1 x syringe A3,5 with 1,0 ml + cannulas
 REF WP4350 1 x syringe B1 with 1,0 ml + cannulas
 REF WP4196 1 x syringe B2 with 1,0 ml + cannulas
 REF WP4398 1 x syringe C2 with 1,0 ml + cannulas
 REF WP4370 1 x syringe A1 with 3,5 g + cannulas
 REF WP4380 1 x syringe A1 with 3,5 g + cannulas
 REF WP4371 1 x syringe A2 with 3,5 g + cannulas
 REF WP4381 1 x syringe A2 with 3,5 g + cannulas
 REF WP4372 1 x syringe A3 with 3,5 g + cannulas
 REF WP4382 1 x syringe A3 with 3,5 g + cannulas
 REF WP4373 1 x syringe A3,5 with 3,5 g + cannulas
 REF WP4383 1 x syringe A3,5 with 3,5 g + cannulas
 REF WP4384 1 x syringe B1 with 3,5 g + cannulas
 REF WP4385 1 x syringe B2 with 3,5 g + cannulas

Disposal

Unused material can be disposed of with household waste after curing.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

We know how

dent a pharm
Produktionsgesellschaft mbH
 Schusterring 35
 25355 Barmstedt/Germany
 Fon: +(49) 4123/9225-0
 Fax: +(49) 4123/9225-48
 www.dent-a-pharm.de
 info@dent-a-pharm.de

Willmann & Pein GmbH
 Schusterring 35
 25355 Barmstedt/Germany
 Fon: +(49) 4123/9228-0
 Fax: +(49) 4123/9228-49
 www.wp-dental.de
 info@wp-dental.de

Made in Germany

H317 May cause allergic skin reactions
 H318 Causes serious eye damage

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