HR

Securafix®

Uputstvo za upotrebu

Securafix je stakloionomerni cement sa posebno dobrim prijanjanjem (adhezijom) na zubni emajl i dentin. Materijal se odlikuje veoma malom topivosti u ustima i maloj kiselosti.

Komponente

Silikatni puder, Poliakrilna kiselina

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Za pričvršćivanje kruna, mostova, inleja, onleja i ortodoksnih traka.

Kontraindikacije

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Securafix®

Instructions for use

Securafix is a glass ionomer cement for attachments with exceptionally good adhesion on tooth enamel and dentine. Its outstanding features include low solubility in the mouth and low acidity.

Components

Silicate powder, Polyacrylic acid

Fields of Application

For attachments of crowns, bridges, inlays, onlays, pins and orthodontic bands.

Contra indications

No direct application on the pulp near dentine.

Method of use

Clean enamel, dentine and metal surfaces carefully and dry thoroughly. It is important to ensure that no traces of eugenol are left which will detract from Securafix effectiveness. It is therefore advisable to attach temporary work with P-Cem eugenol-free cement. Shake powder well before use! Securafix has to be applied at room temperature (15-23°C). In case of higher temperatures is recommended to use a cooled glass-plate. Put powder and liquid on a glass-plate or on a mixing-pad and work the powder portion by portion into the liquid. Use 2 drops of liquid for 1 level measuring scoop of powder. The mixing time is approx. 30 seconds. Securafix can be worked for approx. 1'20" minutes at room temperature and sets in the mouth within approx. 3'00" minutes. Higher temperatures will shorten the working time and lower temperatures prolong it.

In case of bridges and crowns which are to narrow it is important to ensure that the excess of Securafix is not too great, nor its viscosity too high, since glass ionomer cements have specific flow properties.

Securafix has the right consistency if it

can be applied to the tooth stumps and inside the crown by using a stiff brush. During application and setting time prevent contact of work with water and saliva.

Side effects

Side effects with this medical product are extremely rare if it is properly prepared and applied. Should you encounter undesirable side-effects - even in case of doubt- then please contact us and give a precise description of the circumstances and symptoms.

Processing instructions

Ensure that the processing instructions and precautionary measures described in the supplied product information are always strictly observed. Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and give a precise description of the circumstances.

Supplementary indications

Rinse metal instruments with cold water immediately after use and before any cement still attached has time to set. After removal of orthodontic bands, cement residues can be removed gently and rapidly by means of a ultrasonic scaler. Close powder container immediately after use.

Securafix corresponds to the requirements of DIN EN ISO 9917-1:2004-04.

Presentation

WP4110 24g powder, 10g liquid
WP4093 7g powder, 5g liquid
WP4146 10g liquid

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